

# The 500 Watt 20kHz Ultrasonic Process System

## EQUIPMENT SPECIFICATION



**Model L500/5 – 20**

**Sonic Systems Limited**

### **1.0 GENERAL**

The high intensity ultrasonic liquid processing system Model L500/5-20 has been specifically designed for both research and production work in the treatment of liquids and gasses by ultrasound. It has the advantage that acoustic as well as electrical parameters may be monitored thereby giving the user full process information.

Acoustic energy is coupled to the treatment sample by means of titanium horns and can be controlled by an adjustment of electrical power to the transducer and/or by selecting particular horn geometry.

Vibrational amplitude at the working face of the horn, electrical power to the transducer, acoustic power transmitted to the test sample and operational frequencies are all continually monitored.

Transducer resonance is maintained by means of a true motional feedback network and the transducer amplitude is displayed on a digital meter situated on the front panel. The operational frequency is also monitored and displayed.

A wattmeter measures the ultrasonic power to the transducer and it is also possible to measure the acoustic power transmitted to the treatment sample. Power is displayed on a digital meter. The electronic system comprises a main chassis which houses the ultrasonic generator and monitoring facilities.

The system can also be controlled and monitored from a remote station using RS232 serial data communication.

### **2.0 CONTROLS**

- |                              |  |
|------------------------------|--|
| <b>2.1. Mains On/off</b>     | By means of a single-pole rocker switch located on the rear panel. |
| <b>2.2. Sonics On/off</b>    | Switch - Green LED indicates sonics on                             |
| <b>2.3. Amplitude Set</b>    | Variable control on front panel Generator.                         |
| <b>2.4. Local/Remote Set</b> | Switch - LEDs indicate Mode.                                       |



## **6.0 MONITORS**

### **6.1. Amplitude**

Digital meter calibrated in microns peak-peak and measuring transducer displacement. This is presented as a numerical value & bar graph.

Numeric : Microns

Bar : 0-100 %

### **6.2. Frequency**

Digital display measuring operational frequency and Calibrated in Hertz

### **6.3. Power**

Digital meter on front panel measuring Ultrasonic power to the transducer. This is presented as a numerical value & bar graph.

Numeric : Watts

Bar : 0-100 %

## **7.0 ACOUSTIC SYSTEM**

### **7.1. Transducer**

A 20kHz half wavelength Piezo-electric sandwich transducer constructed from Titanium alloy. This is mounted in a fan cooled housing which enables continuous reliable operation.

### **7.2. Frequency**

20kHz nominal

### **7.3. Transducer Displacement**

The maximum vibrational amplitude at the transducer radiating face is 12 +/- 1 microns peak to peak.

### **7.4. Process Horn**

Titanium horns can be supplied for use with each system, giving a variety of displacements and total power output characteristics

## **Sonic Systems Limited**

## **8.0 ACOUSTIC POWER TRANSFER**

There are a range of standard acoustic horn configurations available however coupling high powers into liquid loads can be complex, and an understanding of the acoustic parameters of the load will enable the system performance to be maximised. Sonic Systems can design and manufacture horns and inter-stages for specific application Requirements.

Please refer to our website for reference on considered design :

<http://sonicsystems.co.uk/page/power-ultrasonics-a-guide>

## **9.0 ELECTRICAL**

**9.1. Mains Supply** 230 volts +/-7%, 50 - 60 Hz Single Phase.

**9.2. Classification** Class 1. Must be earthed.

**9.3. Power Consumption** 720VA Maximum

**9.4. Mains connector** IEC Inlet mounted on rear panel.  
Contains integral rocker switch & 20mm fuse.

**9.5. Fuse** 5 Amp. 20mm Ceramic body.

**9.6. Operational Frequency** 20kHz nominal

### **9.7. Frequency Control**

Automatic. A feedback network senses the motion of the transducer and maintains the system at its optimum operating frequency.

## 9.8. Power Output

The generator is designed for a potential maximum output of 500 Watts.

This is a power by demand system where the required transducer displacement can be pre-set and the generator automatically delivers the necessary power to maintain this displacement under varying load conditions.

## 9.9. Duty cycle

Continuous. Capable of 24 hour operation at the approved system power output level

## 10.0 MECHANICAL

The Generators can be easily adapted for 19" rack mounting if required.

### 10.1. Generator Dimensions

Width 444mm  
Depth 350mm  
Height 187mm

### 10.2. Colour

Frame RAL 7030 Grey  
Covers RAL 7035 Light Grey

### 10.3. Transducer Enclosure

Diameter: 65mm (housing)  
Length including connections: approx. 162mm

### 10.4. Weight

Generator Assembly : 13kg  
Transducer Assembly: 925g

## 11.0 DESIGN STANDARDS



### **IEC 1010-1**

Title: Safety requirements for electrical equipment for measurement, control, and laboratory use-Part 1

### **BS.EN 60204-1:1993**

Title: Safety of machinery. Electrical equipment of machines. Specification for general requirements.

### **73/23/EEC Amended 93/68/EEC**

Title: Low Voltage Directive

**12.0 EMC**


Emissions: BSEN61000-6-4 :2001  
Immunity: BSEN61000-6-2: 2005

**13.0 Environmental**

- 13.1. Generator Cabinet** IP20
- 13.2. Transducer assembly** IP20
- 13.3. Treated Liquid Temperature** Max 75 °C
- 13.4. Ambient Temperature** Rated range of use +5°C to+ 35°C
- 13.5. Storage and transport** -40°C to +70 °C

**Please Note:**

The contents of this equipment specification are believed to be correct at the time of printing. The manufacturers, however, reserve the right to change the content, product specification, and performance criteria, without notice.



**Sonic Systems Limited**